

▶ **Naso-Pharyngo-
Laryngoskop 2.8/3.2/3.4**
Naso-Pharyngo-Laryngoscope 2.8/3.2/3.4



grenzenlos im Detail! ▶▶





Naso-Pharyngo-Laryngoskop 2.8/3.2/3.4



► Visionary Medical Technology

Naso-Pharyngo-Laryngoscope 2.8/3.2/3.4



Technische Daten Technical Data

Model	/ Model	FNS-2800	FNS-3200	FNS-3400
Einführschlauch	/ insertion tube	Ø 2,8mm	Ø 3,2mm	Ø 3,4mm
Nutzlänge	/ working length	300mm	300mm	300mm
Blickwinkel	/ field of view	85°	85°	85°
Abwinklung auf/ab angulation up/down		160°/160°	160°/160°	160°/160°
Zubehör	/ accessories	Dichtigkeitstester, Transportbox leakage tester, transportbox	Dichtigkeitstester, Transportbox leakage tester, transportbox	Dichtigkeitstester, Transportbox leakage tester, transportbox
Artikel Nr.	/ item no.	FNS-2800G ■	FNS-3200G ■	FNS-3400G ■
		FNS-2800B ■	FNS-3200B ■	FNS-3400B ■
		FNS-2800R ■	FNS-3200R ■	FNS-3400R ■